



Recombinant Human Serine/arginine repetitive matrix protein 4 (SRRM4)

Product Code	CSB-EP012292HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A7MD48
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MASVQQGEKQ LFEKFWRGTF KAVATPRPES IIVASITARK PLPRTEPQNN PVVPAQDGPS EKLQHLATE PLGTNSWERD KTCRELGATR GHSASHDKDL TPPPSSRGKK KKKKSTRKKR RRSSSYSPSP VKKKKKKSSK KHKRRRSFSK KRRHSSSSPK SKRRDEKRHK KQSRSRPRKS HRHRHRCPS RSQSSSRPS SCESRHRGRS PEEGQKSRRR HSRRCSKTLK KDSPEAQSSR PPSQPLQMLG YLSARGVITG SGSAADLFTK TASPLTTSRG RSQEYDSGND TSSPPSTQTS SARSRGQEKG SPSGGLSKSR ELNSGNTSDS GNSFTTSSPQ NKGAMLENLS PTSRGRESRG FQSPCLECAE VKKSSLVPST ARSSPMKGCS RSSSYASTRS SSHSSRSPNP RASPRYTQSR STSSEKRSYS RSPSYSSKSG KRSPPSRSSR SRRSPSYSRY SPSRERDPKY SEKDSQQRER ERARRRRRSY SPMRKRRRDS PSHLEARRIT SARKRPIPIYY RPSPSSSGSL SSTSSWYSSS SSRASRSYS RSRRSRSRR RSRTRTSSSS SSRSPSPGSR RSRSRSRSR SRSRSQSRYS SSADSYSSTR R
Source	E.coli
Target Names	SRRM4
Protein Names	Recommended name: Serine/arginine repetitive matrix protein 4 Alternative name(s): Medulloblastoma antigen MU-MB-2.76 Neural-specific serine/arginine repetitive splicing factor of 100 kDa Short name= Neural-specific SR-related protein o
Expression Region	1-611
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.