



Recombinant Zea mays Ribosome-inactivating protein

Product Code	CSB-MP333363ZAX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P28522
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	KRIV PKFTEIFPVE DANYPYSAFI ASVRKDVIKH CTDHKGIFQP VLPPEKKVPE LWFYTELKTR TSSITLAIRM DNLYLVGFRT PGGVWWEFGK DGDTHLLGDN PRWLGFGGRY QDLIGNKGLE TVTMGRAEMT RAVNDLAKKK K
Source	Mammalian cell
Protein Names	Recommended name: Ribosome-inactivating protein EC= 3.2.2.22 Alternative name(s): rRNA N-glycosidase Cleaved into the following 2 chains: 1. Ribosome-inactivating protein alpha chain 2. Ribosome-inactivating protein beta chain
Expression Region	17-161
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.