



Recombinant *Saccharomyces cerevisiae* Chorismate synthase (ARO2)

Product Code	CSB-BP333381SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P28777
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	STFGKLF RV TTYGESHCKS VGCIVDGVPP GMSLTEADIQ PQLTRRRPGQ SKLSTPRDEK DRVEIQSGTE FGKTLGTPIA MMIKNE DQRP HDYSDMDKFP RPSHADFTYS EKYGIKASSG GGRASARETI GRVASGAI AE KFLAQNSNVE IVAFVTQIGE IKMNRDSFDP EFQHLLNTIT REKVDSMGPI RCPDASVAGL MVKEIEKYRG NKDSIGGVVT CVVRNLPTGL GEPCFDKLEA MLAHAMLSIP ASKGFEIGSG FQGVSVPGSK HNDPFYFEKE TNRLRKTNN SGGVQGGISN GENIYFSVPF KSVATISQEQ KTATYDGE EG I LA AKGRHDP AVTPRAIPIV EAMTALVLAD ALLIQKARDF SRSVVH
Source	Baculovirus
Target Names	ARO2
Protein Names	Recommended name: Chorismate synthase EC= 4.2.3.5 Alternative name(s): 5-enolpyruvylshikimate-3-phosphate phospholyase
Expression Region	2-376
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.