



Recombinant Drosophila melanogaster Protein scalloped (sd)

Product Code	CSB-MP333498DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P30052
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MKNITSSSTC STGLLQLQNN LSCSELEVAE KTEQQAVGPG TIPSPWTPVN AGPPGALGSA DTNGSMVDSK NLDVGDMSDD EKDLSSADAE GVWSPDIEQS FQEALSIYPP CGRRKIILSD EGKMYGRNEL IARYIKLRTG KTRTRKQVSS HIQVLARRKL REIQAKIKVQ FWQPGLQPST SQDFYDYSIK PFPQPPYPAG KTSTAVSGDE TGIPPSQLPW EGRAIATHKF RLLEFTAFME IQRDEIYHRH LfvQLGGKPS FSDP LLETVD IRQIFDKFPE KSGGLKDLYE KGPQNAFYLV KCWADLNTDL TtGSETGDFY GVTSQYESNE NVVLVCSTIV CSFGKQVVEK VESEYSRLEN NRYVYRIQRS PMCEYMINFI QKLKLNLPERY MMNSVLENFT ILQVMRARET QETLLCIAYV FEVAAQNSGT THHIYRLIKE
Source	Mammalian cell
Target Names	sd
Protein Names	Recommended name: Protein scalloped
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.