



Recombinant Salmonella typhimurium 4-diphosphocytidyl-2-C-methyl-D-erythritol kinase (ispE)

| | |
|--------------------------|---|
| Product Code | CSB-EP333555SXB-B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P30753 |
| Product Type | Recombinant Protein |
| Immunogen Species | Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MMTHWPSPAK LNLFLYITGQ RADGYHTLQT LFQFLDYGDT LHIEPRHDGE IHELLTPVNGV ENEDNLIVRA ARLLMKVASE SGRLPAGSGA DISIEKRLPM GGGLGGSSN AATVLVALNH LWQCGLSIDE LATLGLTLGA DVPV FVRGHA AFAEGVGEIL TPVNPPEKWY LVAHPGVSIP TPVIFKDPQL PRNTPKRSID TLLKCEFSND CEVIARKRFR EVDAALSWLL EYAPSRLTGT GACVFAEFD T ESCARQVLEQ APEWLNFAFVA KGVNLSPLHR ELL |
| Source | E.coli |
| Target Names | ispE |
| Protein Names | Recommended name: 4-diphosphocytidyl-2-C-methyl-D-erythritol kinase Short name= CMK EC= 2.7.1.148 Alternative name(s): 4-(cytidine-5'-diphospho)-2-C-methyl-D-erythritol kinase |
| Expression Region | 1-283 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |