



Recombinant Uncharacterized anhydro-N-acetylmuramic acid kinase-like protein

Product Code	CSB-BP333607YAAQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P31485
Product Type	Recombinant Protein
Immunogen Species	Yersinia enterocolitica
Purity	>85% (SDS-PAGE)
Sequence	MGNKVIAMKS GRFIGVMSGT SLDGIDVVLA AIDERMVAQQ ASYCHPMPLQ LKKDILGMCQ GQSTTL SAVG KLESQLGILF AEAVLALLAK AGLTAQDITA IGCHGQTVWH EPLGEP AFTM QLGDNNRIAA MTKIATVGDF RRRDMAYGGQ GAPLVPAFHH ALLAHSTERR MVLNIGGIAN LSMLLPDLPV RGFDTGPGNM LMDAWIWRNR SLPYDKDNEM DAPTFTASEC LNVSVKTQKK ERPHESIYSW YRLGKKCFPA SWCRLQRSNC S
Source	Baculovirus
Protein Names	Recommended name: Uncharacterized anhydro-N-acetylmuramic acid kinase-like protein
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.