



Recombinant Escherichia coli Thiamine import ATP-binding protein ThiQ (thiQ)

Product Code	CSB-BP333615ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P31548
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MLKLTDITWL YHHLPMRFSL TVERGEQVAI LGPSGAGKST LLNLIAGFLT PASGSLTIDG VDHTTMPPSR RPVSMLFQEN NLFSLTVAQ NIGLGLNPGL KLNAVQQGKM HAIARQMID NLMARLPGEL SGGQRQRVAL ARCLVREQPI LLLDEPFSAL DPALRQEMLT LVSTSCQQQK MLLMVSHSV EDAARIATRS VVVADGRIAW QGMTNELLSG KASASALLGI TG
Source	Baculovirus
Target Names	thiQ
Protein Names	Recommended name: Thiamine import ATP-binding protein ThiQ EC= 3.6.3.-
Expression Region	1-232
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.