



Recombinant Drosophila melanogaster Fork head domain transcription factor slp1 (slp1)

Product Code	CSB-MP333661DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32030
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MVKSEMEFKS NFSIDAILAK KPINTATQPI KTEPVHHHHQ YVHPYSNSDG ELSASEDFDS PSRTSTPMSS AAESLSSQNN DKLDVEFDDE LEDQLDEDEQE SEDCGNPSKKQ KMTAGSDTKK PPYSYNALIM MAIQDSPEQR LTLNGIYQYL INRFPYFKAN KRGWQNSIRH NLSLNKCFTK IPRSYYDDPGK GNYWILDPSA EEVFIGETTQ KLRRKNPGAS RTRLAAYRQA IFSPMMAASP YGAPAPSYGY PAVPFAAAAA AALYQRMNPA AYQAAYQQMQ YQQAPQAHHH QAPHPAQMQG YPPQLNAELF QRMQFFGKFP SS
Source	Mammalian cell
Target Names	slp1
Protein Names	Recommended name: Fork head domain transcription factor slp1 Alternative name(s): Sloppy paired locus protein 1
Expression Region	1-322
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.