



Recombinant Escherichia coli Putative uncharacterized protein yehQ (yehQ)

Product Code	CSB-MP333788ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P33353
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MNSLRPELLE LTPQALTALS NAGFVKRSLK ELENGNVPEI SHENDALIAT FSDGVRTQLA NGQALKEAQC SCGANGMCRH RVMLVLSYQR LCATTQSTEK EEEWDPAIWL EELATLPDAT RKRAQALVAK GITIELFCAP GEIPSARLPM SDVRFYSRSS IRFARDCIE GTLCEHVLA VQAFVEAKAQ QAEFNHLIWQ MRSEHVTSSD DPFASEEGNA CRQYVQQLSQ TLWLGGISQP LIHYEAAFNR ALQAAETCNW RWVSESLRQL RASVDAFHAR ASHYNAGECL HQLAALNSRL NCAQEMARRD SIGEVPPVPW RTVVGSGIAG EAKLDHLRLV SLGMRCWQDI EHYGLRIWFT DPDTGSILHL SRSWPRSEQE NSPAATRRLF SFQAGALAGG QIVSQAARKS ADGELLATR NRLSSVPLS PDAWQMLSAP LRQPGIVALR EYLRQRPPAC IRPLNQVDNL FILPVAECIS LGWDSSRQTL DAQVISGEGE DNVLTLSLPA SASAPYAVER MAALLQQTDD PVCLVSGFVS FVEGQLTLEP RVMMTKTRAW ALDAETTPVA PLPSASVLPV PSTAHQLLIR CQALLIQLLH NGWRYQEQA IGQA
Source	Mammalian cell
Target Names	yehQ
Protein Names	Recommended name: Putative uncharacterized protein yehQ
Expression Region	1-614
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.