



Recombinant Bovine Transcription factor jun-D (JUND)

Product Code	CSB-MP011974BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A7YY54
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	METPFYGDEA LSGLGGGGSS SGGGGSFASP GRLFPGAPPT AAPGSMKKD ALTLSEQV AAALKPAAAP PPGPLRTDGA PGTAPPDGLL ASPELGLLKL ASPELERLII QSNGLVTTTP TSTQFLYPKV AASEEQEFAE GFVKALEDLH KQNQLSAGAA SAAAAAGGPS GTAAGAAPPS ELAPAAATPE APVYANLSSY AGGTGSAGGA ATVAFAAEPV PFPPPPPGT LGPPRLAALK DEPQTPDVP SFGESPPLSP IDMDTQERIK AERKRLRNRI AASKCRKRKL ERISRLEEKV KTLKSQNTL ASTASLLREQ VAQLKQKVLV HVNSGCQLLP QHQPAY
Source	Mammalian cell
Target Names	JUND
Protein Names	Recommended name: Transcription factor jun-D
Expression Region	1-347
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	The protein encoded by this intronless gene is a member of the JUN family, and a functional component of the AP1 transcription factor complex. It has been proposed to protect cells from p53-dependent senescence and apoptosis. Alternate translation initiation site usage results in the production of different isoforms.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.