



Recombinant ATP-dependent Clp protease ATP-binding subunit ClpX (clpX)

Product Code	CSB-EP333816DPX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P33683
Product Type	Recombinant Protein
Immunogen Species	Azotobacter vinelandii
Purity	>85% (SDS-PAGE)
Sequence	MDKTDKPCS FCGAEKSPTV PLIAGNDGRI CEACVKLAHQ VVSSWGQRRK NQQLAPQLRT PVEYKKHLDE SVIGQEAKE TLSVAVYNHY LRLNCDREP VCQLGGEQVEL EKSNI LMAGP SGTGKTL LVR TLARILGV PF AMADATTLTQ AGYVGDDVDS IITRLLDAAG GDVQQAQWGI VYIDEVDKLA KRSGGGTAVR DISGEGVQQA LLKMVEGTEV RISKSGRRNE HSEEQVVDTR NILFIAGGAF PGLEALVSSR IQPKNTGIGF HAQPRREAPS INELMASLLP DDLHEFGLIP EFIGRFPIIT FLQELDHATL LRILTEPRNA LVKQYKQLFA YQGVELVITD AALNYIADQA LIRKTGARGL RAVLEAALQQ TMFNMPSPQP LRGCTLDLVE HEDGGRSLEV LTRLAEDGSG RIQPDPSPVV EEKSALSADL
Source	E.coli
Target Names	clpX
Protein Names	Recommended name: ATP-dependent Clp protease ATP-binding subunit ClpX
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.