



Recombinant Arabidopsis thaliana Proteasome subunit alpha type-1-A (PAF1)

Product Code	CSB-EP333854DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34066
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MFRNQYDTDV TTWSPTGRLF QVEYAMEAVK QGSAAIGLRS RSHVVLACVN KAQSELSSHQ RKIFKVDHDI GVAIAGLTAD GRVLSRYMRS ESINHSFTYE SPLPVGRLVV HLADKAQVCT QRSWKRPYGV GLLVGGLDES GAHLYYNCPS GNYFEYQAF AIGSRSQAAKT YLERRFESFG DSSREDLIKD AILAVRETLQ GETLKSSLCT VAILGVDEPF HFLDQEAIQK VIDTFEKVPE EEEGEGEAGE GEAEAAEAAP AERGGGVAGD QDVAPMEM
Source	E.coli
Target Names	PAF1
Protein Names	Recommended name: Proteasome subunit alpha type-1-A EC= 3.4.25.1 Alternative name(s): 20S proteasome alpha subunit F-1 Proteasome 30 kDa subunit Proteasome component 2A Short name= AtPSM30 Proteasome subunit alpha type-6
Expression Region	1-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.