



Recombinant Escherichia coli Outer membrane protease ompP (ompP)

Product Code	CSB-MP333866ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34210
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	SDFFGPE KISTEINLGT LSGKTKERVY EPEEGGRKVS QLDWKYSNAA ILKGAVNWEL NPWLSVGAAG WTTLNSRGGN MVDQDWMDSG TPGTWTDESR HPDTRLNYAN EFDLNVKGMWF LKESDYRLAI MAGYQESRYS FNATGGTYIY SENGGFRNET GALPDKIKVI GYKQHFPIPY VGLTGNYRYD NFEFGGAFKY SGWVRGSDND EHYVRQTTFR SKVINQNYYS VAVNAGYYIT PEAKVYIEGV WSRLTNKKGD TSLYDRSDNT SEHNNNGAGI ENYNFITTAG LKYTF
Source	Mammalian cell
Target Names	ompP
Protein Names	Recommended name: Outer membrane protease ompP EC= 3.4.23.-
Expression Region	24-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.