



Recombinant Uncharacterized F-box/LRR-repeat protein C02F5.7 (C02F5.7)

Product Code	CSB-EP333876CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34284
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	METAERISAA ASAASSRRAK RLAQQAHKTH PVIQAKQNQM YLITTLSPAQ VDNSLINRVL PKEVLLKVFS FLDTKALCRS AQVCRSWSIL ALDGSNWQRV DLFTFQRDVK TAVVENLARR CGGFLKELSL KGCENVHDSA LRTFTSRCPP LEHLSLYRCK RVTDASCENL GRYCHKLNLYL NLENCSSITD RAMKYIGDGC PNLSYLNISW CDAIQDRGVQ IILSNCKSLD TLILRGCEGL TENVFGSVEA HMGAIKKLNL LQCFQLTDIT VQNIANGATA LEYLCMSNCN QISDRSLVSL GQSHNLKVL ELSGCTLLGD NGFIPLARGC RQLERLDMED CSLISDHTIN SLANNCTALR ELSLSHCELI TDESIQNLAS KHRETLNVLE LDNCPQLTDS TLSHLRHCKA LKRIDLDCQ NVSKEAIVRF QHHRPNIEIH AYFAPVTPPT DQVVNRGGIC RCCVIL
Source	E.coli
Target Names	C02F5.7
Protein Names	Recommended name: Uncharacterized F-box/LRR-repeat protein C02F5.7
Expression Region	1-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.