



Recombinant Uncharacterized F-box protein C05B5.5 (C05B5.5)

Product Code	CSB-YP333877CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34293
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MGQYLTKTEE SSLNYVDLPD TVHRKIFEYL NPWEIFKLSR ISKAIHVTIL KNKKFAVKDI DLCTDENILK FQFQFVNNIV LSWEFYNLHE KPYFTSQFSN RCQYRQQYDI KNYENPEEA FLFAFRNALS IFNVQNSKLN RFYLAPSKLE LFFKLPYSFI DGRNCENLSI WNKKANDEDE RDFSINFAKA VLEKMGVVRN LRLLFNCSTK DNEFDCEKEY YNSETDKLMP NDWYLYNKFD DIDLNVFLKN WLNFRDPLIR KFRIRGKRNF NTAIFDGIE TIPWDQRNGL NVYYYSATAN DTGSLIFQAN ENHAIVKVGGA INGNESGYM EFLFITF
Source	Yeast
Target Names	C05B5.5
Protein Names	Recommended name: Uncharacterized F-box protein C05B5.5
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.