



# Recombinant Bovine Keratin, type I cytoskeletal 40 (KRT40)

<b>Product Code</b>	CSB-BP012557BO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A7YWM2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MASDCSPTGC SSESSARASD CALASTCSVE TTCLPSACAT SSCQTPSFPS GARLPTGCLP PACFAGSCNI PCVVGNCWC EDGVFNSNEK ETMQFLNDR ASYLEKVRGL EELNAELECR IREQCEEDVP LVCPDYQCYF DTIEDLQQKI LCTKAENCRL AVQLDNCKLA ADDFRSKYES ELSLRQLVET DISGLRGILG ELTVCRSDLE AHVESLKDDL LCLKSHEEE VNVLRGQLGD RLSVELDTAP TTDLNRVLDE MRCQYETVLA NNRRDVEEWF AAQTEELNQQ QLSSAEQLQG CQTEILELKR TANTLEIELQ AQQSLTESLE CTVAETEAQY SSELAQIQCL IDNVENQLAE IRCDLERQNG EYRVLLDTKA RLECEINTYQ GLLDSEDSRL PCNPCSATSM SNTDCEPCSA YVICTVENSC P
<b>Source</b>	Baculovirus
<b>Target Names</b>	KRT40
<b>Protein Names</b>	Recommended name: Keratin, type I cytoskeletal 40 Alternative name(s): Cytokeratin-40 Short name= CK-40 Keratin-40 Short name= K40
<b>Expression Region</b>	1-431
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This gene encodes a member of the type I (acidic) keratin family, which belongs to the superfamily of intermediate filament (IF) proteins. Keratins are heteropolymeric structural proteins which form the intermediate filament. These filaments, along with actin microfilaments and microtubules, compose the cytoskeleton of epithelial cells. The type I keratin genes are clustered in a region of chromosome 17q12-q21.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.