



Recombinant Human adenovirus B serotype 11 Fiber protein (L5)

Product Code	CSB-MP334023HID
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P35773
Product Type	Recombinant Protein
Immunogen Species	Human adenovirus B serotype 11 (strain BC34) (HAdV-11) (Human adenovirus 11A (strain BC34))
Purity	>85% (SDS-PAGE)
Sequence	MTKRVRLSDS FNPVYPYEDE STSQHPFINP GFISPNGFTQ SPDGVLTLC LTPLTTTGGG LQLKVGGGLT VDDTDGTLQE NIGTTTPLVK TGHSIGLSLG AGLGTDENKL CTKLGKGLTF NSNNICIDN INTLWTGINP TEANCQMMDS SESNCKLIL TLVKTGALVT AFVYVIGVSN NFNMLTTYRN INFTAELFFD SAGNLLTSL SLKTPLNHKS GQNMATGAI NAKSFMPSTT AYPFNNNSRE KENYIYGTCH YTASDHAFP IDISVMLNQR AIRADTSYCI RITWSWNTGD APEGQTSATT LVTSPFTFY IREDD
Source	Mammalian cell
Target Names	L5
Protein Names	Recommended name: Fiber protein Alternative name(s): pIV
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.