



Recombinant *Saccharomyces cerevisiae* Tyrosine-tRNA ligase, cytoplasmic (TYS1)

Product Code	CSB-MP334127SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P36421
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	SSAATVDPN EAFGLITKNL QEVLNPQIIK DVLEVQKRHL KLYWGTAPTG RPHCGYFVPM TKLADFLKAG CEVTVLLADL HAFLDNMKAP LEVVNYRAKY YELTIKAILR SINVPICKLK FVVGSSYQLT PDYTMDIFRL SNIVSQNDK RAGADVVKQV ANPLLSGLIY PLMQALDEQF LDVDCQFGGV DQRKIFVLAE ENLPSLGYKK RAHLMNPMVP GLAQGGKMSA SDPNSKIDLL EEPKQVKKKI NSAFCSPGNV EENGLLSFVQ YVIAPIQELK FGTNHFEFFI DRPEKFGGPI TYKSFEEMKL AFKEEKLSPD DLKIGVADAI NELLEPIRQE FANNKEFQEA SEKGYPVATP QKSKKAKKPK NKGTKYPGAT KTNEIATKLE ETKL
Source	Mammalian cell
Target Names	TYS1
Protein Names	Recommended name: Tyrosine--tRNA ligase, cytoplasmic EC= 6.1.1.1 Alternative name(s): Tyrosyl-tRNA synthetase Short name= TyrRS
Expression Region	2-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.