



Recombinant Escherichia coli Probable endopeptidase YaeF (yaeF)

Product Code	CSB-EP334218ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P37056
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	CTVDISQPDP SATAVDAAEK TWAVKFQHQSFTEQSIKEI TAPDLKPGDL LFSSSLGVTS FGIRVFSTSS VSHVAIFLGD NNVAEATGAG VQIVSLKKAM KHSDKLFVLR VPDLTQQAT DITAFANKIK DSGYNYRGIV EFIPFMVTRQ MCSLNPFSERFRQQCVSGLA KAQLSSVGEG DKKSWFCSEF VTDAFAKAGH PLTLAQSGWI SPADLMHMRI GDVSAFKPET QLQYVGHLKP GIYIKAGRFV GLTR
Source	E.coli
Target Names	yaeF
Protein Names	Recommended name: Probable endopeptidase YaeF EC= 3.4.-.- Alternative name(s): Uncharacterized lipoprotein YaeF
Expression Region	21-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.