



Recombinant Sorbose-specific phosphotransferase enzyme IIB component (sorB)

Product Code	CSB-YP334223KBG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P37081
Product Type	Recombinant Protein
Immunogen Species	Klebsiella pneumoniae
Purity	>85% (SDS-PAGE)
Sequence	MQITLARIDD RLIHGQVTTV WSKVANAQRI IICNDDVFND EVRRTLLRQA APPGMKVNVV SLEKAVAVYH NPQYQDETVF YLFTNPHDVL TMVRQGVQIA TLNIGGMAWR PGKKQLTKAV SLDPQDIQAF RELDKLGVKL DLRVVASDPS VNILDKINET AFCE
Source	Yeast
Target Names	sorB
Protein Names	Recommended name: Sorbose-specific phosphotransferase enzyme IIB component EC= 2.7.1.69 Alternative name(s): EIIB-Sor EIII-B-Sor PTS system sorbose-specific EIIB component
Expression Region	1-164
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.