



Recombinant Escherichia coli HTH-type transcriptional repressor CsiR (csiR)

Product Code	CSB-YP334250ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P37338
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MTITSLDGYR WLKNDIIRGN FQPDEKLRMS LLTSRYALGV GPLREALSQL VAERLVTVVN QKGYRVASMS EQELLDIFDA RANMEAMLVS LAIARGGDEW EADVLAKAHL LSKLEACDAS EKMLDEWDLR HQAFHTAIVA GCGSHYLLQM RERLFDLAAR YRFIWLRRTV LSVEMLEDKH DQHQLTAAV LARDTARASE LMRQHLLTPI PIIQQAMAGN
Source	Yeast
Target Names	csiR
Protein Names	Recommended name: HTH-type transcriptional repressor CsiR Alternative name(s): Carbon starvation induced regulator
Expression Region	1-220
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.