



Recombinant Escherichia coli Uncharacterized protein yhiJ (yhiJ)

Product Code	CSB-MP334295ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P37627
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MKIGTVAGTN DSTTTIATND MVQEHVTNFT KELFGYIANG IGDDISSIAR TMLGEVVEKI DDWQUIERFQQ SIQDDKISFT IQTDHSEKYS MLSGMRAHIL RRNNNYQFIV TINSKNYGCS LDNTDINWCS IVYLLNNMTV NDNANDVAVT ESYKPIWNWK ISQYNVSDIK FETMIKPQFA DRIYFSNCLP VDPTSTRPTY FGDTDGSVGGA VLFALFATGH LGIMAEGENF LSQLLNIEDE VLNVLLENF NEQLNTNVNT IISILNRRDI ILESLOPYLV INKDAVTPCT FLGDQTGDRF SNICGDQFII DLLKRIMSIN ENVHVLAGNH ETNCNGNYMQ NFTRMKPLDE DTYSGIKDYP VCFYDPKYKI MANHHGITFD DQRKRYIIGP ITVSIDEMTN ALDPVELAAI INKKHHAIIN GKKFKTSRAI SCRSFNRYFS VSTDYRPKLE ALLACSQMLG INQVVAHNGN GGRERIGETG TVLGLNARDS KHAGRMFSMH NCQINPGAGP EITTPWKSQYQ HEKNRNLMP LIRRTMLQL
Source	Mammalian cell
Target Names	yhiJ
Protein Names	Recommended name: Uncharacterized protein yhiJ
Expression Region	1-540
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.