



Recombinant Human WAS protein family homolog 1 (WASH1)

Product Code	CSB-BP025972HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A8K0Z3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTPVVMQHSL AGQTYAVPFI QPDLRREEAV QQMADALQYL QKVSGDIFSR ISQQVEQSRS QVQAIGEKVS LAQAKIEKIK GSKKAIKVFS SAKYPAPGRL QEYGSIFTGA QDPGLQRRPR HRIQSKHRPL DERALQEKLK DFPVCVSTKP EPEDDAEEGL GGLPSNISSV SLLLLFNTTE NLYKKYVFLD PLAGAVTKTH VMLGAETEEK LFDAPLSISK REQLEQQVPE NYFYVPDLGQ VPEIHVPSYL PDLPGIANDL MYSADLGPGI APSAPGTIPE LPTFHTEVAE PLKVDLQDGV LTPPPPPPPP PPAPEVLASA PPLPPSTAAP VGQGARQDDS SSSASPSVQG APREVDPSPG GWATLLESIR QAGGIGKAKL RSMKEREKLEK QQQKEQEQVR ATSQGGHLMs DLFNKLVMRK KGISGKGPGA GEGPGGAFVR VSDSIPPLPP PQQPQAEDEE DDWES
Source	Baculovirus
Target Names	WASHC1
Protein Names	Recommended name: WAS protein family homolog 1 Alternative name(s): CXYorf1-like protein on chromosome 9 Protein FAM39E
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.