



Recombinant Human Methyltransferase-like protein 12, mitochondrial (METTL12)

Product Code	CSB-BP013727HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A8MUP2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	CL ADRCLWDR LH AQPRLGTVPT FDWFFGYDEV QGLLLPLLQE AQAASPLRVL DVGCGTSSLC TGLYTKSPHP VDLGVDFSP VAVAHMNSLL EGGPGQTPLC PGHPASSLHF MHADAQNLGA VASSGSFQLL LDKGTWDAVA RGGLPRAYQL LSECLRVLNP QGTLIQFSDE DPDVRLPCLE QGSYGWTVTV QELGPFRGIT YFAYLIQGS
Source	Baculovirus
Target Names	CSKMT
Protein Names	Recommended name: Methyltransferase-like protein 12, mitochondrial EC=2.1.1.-
Expression Region	29-240
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.