



Recombinant Human E3 ubiquitin-protein ligase NEURL1B (NEURL1B)

Product Code	CSB-BP015722HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A8MQ27
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MGNTVHRTL P DPSPARLLA TRPCCGPGPE RRPVLGEAPR FHAQAKGKNV RLDGHSRRAT RRNSFCNGVT FTQRPIRLYE QVRLRLVAVR PGWSGALRFG FTAHDPSLMS AQDIPKYACP DLVTRPGYWA KALPENLALR DTVLAYWADR HGRVYFSVND GEPVLFHCGV AVGGPLWALI DVYGITDEVQ LLESFADTL TPARLSQARF SACLPPSSHD AANFDNNELE NNQVVAKLGH LALGRAPGPP PADAAAAAIP CGPRERPRPA SSPALLEADL RFHATRGPDV SLSADRKVC APRPDGGRTL VFSERPLRPG ESLFVEVGRP GLAAPGALAF GITSCDPGVL RPNELPADPD ALLDRKEYWV VARAGPVPSG GDALSFTLRP GGDVLLGING RPRGRLLCVD TTQALWAFFA VRGGVAGQLR LLGTLQSSPA TTTPSGSLSG SQDDSDSDMT FSVNQSSAS ESSLVTAPSS PLSPPVSPVF SPPEPAGIKN GECTVCFDGE VDTVYITCGH MCLCHSCGLR LKRQARACCP ICRRPIKDV KIYRP</p>
Source	Baculovirus
Target Names	NEURL1B
Protein Names	Recommended name: E3 ubiquitin-protein ligase NEURL1B EC= 6.3.2.- Alternative name(s): Neuralized-like protein 1B Neuralized-like protein 3
Expression Region	1-555
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.