



Recombinant Escherichia coli Isoaspartyl dipeptidase (iadA)

Product Code	CSB-BP334552ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39377
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MIDYTAAGFT LLQGAHLYAP EDRGICDVLV ANGKIIAVAS NIPSDIVPNC TVVDLSGQIL CPGFIDQHVH LIGGGGEAGP TTRTPEVALS RLTEAGVTSV VGLLGTDSIS RHPESELLAKT RALNEEGISA WMLTGAYHVP SRTITGSVEK DVAIIDRVIG VKCAISDHRS AAPDVYHLAN MAAESRVGGL LGGKPGVTVF HMGDSKKALQ PIYDLENCD VPISKLLPTH VNRNVPLFEQ ALEFARKGGT IDITSSIDEP VAPAEGIARA VQAGIPLARV TLSSDGNQSQ PFFDDEGNLT HIGVAGFETL LETVQVLVKD YDFSISDALR PLTSSVAGFL NLTGKGEILP GNDADLLVMT PELRIEQVYA RGKLMVKDGK ACVKGTFETA
Source	Baculovirus
Target Names	iadA
Protein Names	Recommended name: Isoaspartyl dipeptidase EC= 3.4.19.-
Expression Region	1-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.