



Recombinant *Synechococcus elongatus* Sensor protein sphS (sphS)

Product Code	CSB-EP334598FPY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39664
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	≥85% (SDS-PAGE)
Sequence	MAAWEFALGL LTASLWRWAR KWRSPVKVKP MAAVSSLEP QLEQITDLR DRDRILLEDLP VSFLLLDADN LVLEANRSAR VLLALPPEDY CRPLLEVRS YELDRLVARC RAANAPQTDR WTLTPVNPDP LQVVPQTPRP VQGQAIPLSN GQIGVLIEDR QELVDLAQQR NRWSDVAHE LKTPLOTSIRL LAEALRDRLQ DEPQVWVDRL LGETQRLGQL VQDLLELSRL EQGPSGLQKL EAVDLVALLT SVRNSLEPLA EPLRLGWAYQ GPEQGFVRGD RQRLFRLWLN LVDNAIRHSP SGGCLYVELR QRGDTWICDL YDDGPGFADA DLPYLFERFY RGDPSRVRPA AASSSSPGSG LGLAIARQVV EAHQGRISAR NHPVTGGAWL RVQLPQEPSL TPALKIGTGR RSG
Source	E.coli
Target Names	sphS
Protein Names	Recommended name: Sensor protein sphS EC= 2.7.13.3
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.