



# Recombinant *Saccharomyces cerevisiae* Reduced viability upon starvation protein 167 (RVS167)

<b>Product Code</b>	CSB-MP334610SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P39743
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SFKGFTKAV SRAPQSFRQK FKMGEQTEDP VYEDAERRFQ ELEQETKCLS EESKRYSTAV NGMLTHQIGF AKSMEEIFKP ISGKMSDPNA TIPEDNPQGI EASEQYRAIV AELQETLKPD LALVEEKIVT PCQELLKIIT YIRKMATKRN HKKLDLDRHL NTYNKHEKKEK EPTAKDEERL YKAQAQVEVA QQEYDYYNDL LKTQLPILFS LEAEFVKPLF VSFYFMQLNI FYTLYNRLQD MKIPYFDLNS DIVESYIAKK GNVEEQTDAL TITHFKLGYS KAKLEMTRRK YGVATAEGSP VSGASSGVGY GAGYDPATAT SPTPTGYGYG AAAPSYAAQP AAQYGTAAAV GTAAAVGTAA GAAAGAVPGT YPQYAAAQSP PLTGLGFQQS PQQQQGPPPA YSNPLTSPVA GTPAAAVAAA PGVETVTALY DYQAQAAGDL SFPAGAVIEI VQRTPDVNEW WTGRYNGQQG VFPGNVYQLN KN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	RVS167
<b>Protein Names</b>	Recommended name: Reduced viability upon starvation protein 167
<b>Expression Region</b>	2-482
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.