



Recombinant *Saccharomyces cerevisiae* Protein BIM1 (BIM1)

Product Code	CSB-MP334656SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40013
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	SAGIGESRT ELLTWLNGLL NLNYKKIEEC GTGAAYCQIM DSIYGDLP MN RVKFNATAEY EFQTNKILQ SCFSRHGIEK TVYVDKLIRC KFQDNLEFLQ WLKKHWIRHK DESVYDPDAR RKYRPIITNN SATKPRTVSN PTTAKRSSST GTGSAMSGGL ATRHSSLGIN GSRKTSVTQG QLVAIQAELT KSQETIGSLN EEIEQYKGTV STLEIEREFY FNKL RDIEIL VHTTQDLIN GVKFNDETI TGHGNGNGGA LLRFVKKVES ILYATAEGFE MNDGEDELND KNLGEHGTVP NQGGYANSNG EVNGNEGSNH DVIMQNDEGE VGVSNLIID EETF
Source	Mammalian cell
Target Names	BIM1
Protein Names	Recommended name: Protein BIM1
Expression Region	2-344
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.