



Recombinant Human Guanine nucleotide-binding protein G (t) subunit alpha-3 (GNAT3)

Product Code	CSB-YP009600HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A8MTJ3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	GSGISSESK ESAKRSKELE KKLQEDAERD ARTVKLLLLG AGESGKSTIV KQMKIIHKNG YSEQECMEFK AVIYSNTLQS ILAIVKAMTT LGIDYVNPRS AEDQRQLYAM ANTLEDGGMT PQLAEVIKRL WRDPGIQACF ERASEYQLND SAAYYLNDLD RITASGYVPN EQDVLHSRVK TTGIITQFS FKDLHFRMFD VGGQRSERKK WIHCFEGVTC IIFCAALSAY DMVLVEDEEV NRMHESLHLF NSICNHKYFS TTSIVLFLNK KDIFQEKVTK VHLSICFPEY TGPNTFEDAG NYIKNQFLDL NLKKEDKEIY SHMTCATDTQ NVKVFVDAVT DIIKENLKD CGLF
Source	Yeast
Target Names	GNAT3
Protein Names	Recommended name: Guanine nucleotide-binding protein G(t) subunit alpha-3 Alternative name(s): Gustducin alpha-3 chain
Expression Region	2-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.