



Recombinant Human Serpin E3 (SERPINE3)

Product Code	CSB-MP021083HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A8MV23
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	HLREGMTLLK TEFALHLYQS VAACRNETNF VISPAGVSLP LEILQFGAEG STGQQLADAL GYTVHDKRVK DFLHAVYATL PTSSQGTEME LACSLFVQVG TPLSPCFVEH VSWWANSSE PADLSEPNST AIQTSEGASR ETAGGGPSEG PGGWPWEQVS AAFAQLVLVS TMSFQGTWRK RFSSTDTQIL PFTCAYGLVL QVPMMHQTTE VNYGQFQDTA GHQVGVLELP YLGSVAVSLFL VLPRDKDTPL SHIEPHLTAS TIHLWTTSLR RARMDVFLPR FRIQNQFNLK SILNSWGVTD LFDPLKANLK GISGQDGFYV SEAIHKAKIE VLEEGTKASG ATALLLKRS RIPIFKADRP FIYFLREPNT GITVFFDRIQ IYQCLSSNK GSFVHYPLKN KHSF
Source	Mammalian cell
Target Names	SERPINE3
Protein Names	Recommended name: Serpin E3
Expression Region	21-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.