



Recombinant *Saccharomyces cerevisiae* Nascent polypeptide-associated complex subunit beta-2 (BTT1)

Product Code	CSB-BP334693SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40314
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MPVDQEKLAK LHKLSAANKV GGTRRKINKK GNLYNNNDKD NTKLQAEHLK LHPMTIENVA EANFFKNGK VLHFNSAVVQ IAPQCNTMI HGQPKENTLN GLYPSVASQL GSQELEYLTG LAHNLENEQT VLDQLGDRCS ETKQQVMNS
Source	Baculovirus
Target Names	BTT1
Protein Names	Recommended name: Nascent polypeptide-associated complex subunit beta-2 Short name= NAC-beta-2 Alternative name(s): BTF3 homolog BTT1 Beta-2-NAC
Expression Region	1-149
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.