



Recombinant *Saccharomyces cerevisiae* Autophagy-related protein 3 (ATG3)

Product Code	CSB-BP334699SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40344
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MIRSTLSSWR EYLTPITHKS TFLTGGQITP EEFVQAGDYL CHMFPTWKWN EESDISYRD FLPKNKQLI IRKVPCDKRA EQCVEVEGPD VIMKGFAEDG DEDDVLEYIG SETEHVQSTP AGGTDSSID DIDELIQDME IKEEDENDDT EEFNAKGGLA KDMAQERYD LYIAYSTSYR VPKMYIVGFN SNGSPLSPEQ MFEDISADYR TKTATIEKLP FYKNSVLSVS IHPCKHANVM KILLDKVRVV RQRRRKELQE EQELDGVGDW EDLQDDIDDS LRVDQYLIVF LKFITSVTPS IQHDYTMEGW
Source	Baculovirus
Target Names	ATG3
Protein Names	Recommended name: Autophagy-related protein 3 Alternative name(s): Autophagy-related E2-like conjugation enzyme ATG3
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.