



Recombinant *Saccharomyces cerevisiae* Histidine protein methyltransferase 1 (HPM1)

Product Code	CSB-BP334722SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40481
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSFSFGFTSN DFDDDELVAQ PETFVESSKE NENTTAYINP LDSDFLSQAG VVQPNVEDLG TILES�KDVR LTFEEFQSPI YRKPLIKREL FDVKHQLMLE TDAQSNNNST ELDILLGDTSEDLRKNIEYEG GLKSWECYSYD LVDLLSENVD RISNDIDAVV EIGCGTALPS EFLFRSALLR NDRSKGLKFV LTDYNASVLR LVTIPNLVIT WAKTVLTKEQ WYALQKDECE DIPINNEELL LTKLLAAFY DDVQSRNISV TLISGSWGRK FSNLIHEVLS GSQKVLSSLSS ETIYQPDNLP VIAETILDIH NLPQTDVKTY VAAKDIYFGV GGSITEFEAY LDDKINSEHL PIHSERFKVN SGLKRSIICI ETNKAIR
Source	Baculovirus
Target Names	HPM1
Protein Names	Recommended name: Histidine protein methyltransferase 1 EC= 2.1.1.85 Alternative name(s): Mitotic exit network interactor 1
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.