



# Recombinant *Saccharomyces cerevisiae* DBF2 kinase activator protein MOB1 (MOB1)

<b>Product Code</b>	CSB-YP334723SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P40484
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSFLQNFHIS PGQTIRSTRG FKWNTANAAN NAGSVSPTKA TPHNNTINGN NNNANTINNR ADFTNPNVNG YNESDHGRMS PVLTTPKRHA PPPEQLQNVT DFNYTPSHQK PFLQPQAGTT VTTHQDIKQI VEMTLGSEGV LNQAVKLPRG EDENEWLAVH CVDFYNQINM LYGSITEFCS PQTCPRIAT NEYEYLWAFQ KGQPPVSVSA PKYVECLMRW CQDQFDDES LFPSKVTGTFP EGFIQRVIQP ILRRLFRVYA HIYCHHFNEI LELNLQTVLN TFRHFCLFA QEFELLRPAD FGPLLELVME LRDR
<b>Source</b>	Yeast
<b>Target Names</b>	MOB1
<b>Protein Names</b>	Recommended name: DBF2 kinase activator protein MOB1 Alternative name(s): MPS1 binder 1 Maintenance of ploidy protein MOB1
<b>Expression Region</b>	1-314
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.