



Recombinant Saccharomyces cerevisiae U2 small nuclear ribonucleoprotein B' (MSL1)

Product Code	CSB-EP334740SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40567
Product Type	Recombinant Protein
Immunogen Species	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MVEPARKKQR IDRDTHHTVA EPVTEAKNTL YVSQLNEKIN MQLLRVNLFL LFATFGEVLK VSMNFKKQRG QAFITMRTID QASLAQISLN GERFFGKPLK VEFSKSETKT L
Source	E.coli
Target Names	MSL1
Protein Names	Recommended name: U2 small nuclear ribonucleoprotein B" Short name= U2 snRNP B"
Expression Region	1-111
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.