



Recombinant Bacillus subtilis Aryl-phospho-beta-D-glucosidase BglH (bglH)

Product Code	CSB-MP334762BRJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40740
Product Type	Recombinant Protein
Immunogen Species	Bacillus subtilis (strain 168)
Purity	>85% (SDS-PAGE)
Sequence	MSSNEKRFPE GFLWGGAVAA NQVEGAYNEG GKGLSTADVS PNGIMSPFDE SMTSLNLYHN GIDFYHRYKE DIALFAEMGF KAFRTSIAWT RIFPNGDEEE PNEEGLRFYD DLFDELLKHH IEPVVTISHY EMPLGLVKNY GGWKNRKYVIE FYERYAKTVF KRYQHKVKYW MTFNEINVVL HAPFTGGGLV FEEGENKLNA MYQAAHHQFV ASALAVKAGH DIIPDSKIGC MIAATTTYPM TSKPEDVFAA MENERKTLFF SDVQARGAYP GYMCRYLAEN NIEIEMAEGD EELLKEHTVD YIGFSYYMSM AASTDPEELA KSGGNLLGGV KNPYLKSSEW GWQIDPKGLR ITLNTLYDRY QKPLFIVENG LGAVDKVEED GTIQDDYRIN YLRDHLIAR EAIADGVELI GYTSWGPIDL VSASTAEMKK RYGFIIYVDRD NEGNLTFNRI KKKSFNWYQQ VIATNGESL
Source	Mammalian cell
Target Names	bglH
Protein Names	Recommended name: Aryl-phospho-beta-D-glucosidase BglH EC= 3.2.1.86 Alternative name(s): 6-phospho-beta-glucosidase
Expression Region	1-469
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.