



Recombinant Human RING finger protein C14orf164 (C14orf164)

Product Code	CSB-MP003163HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A8MTL3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDWFHCNQCF RKDGAHFFVT SCGHIFCKKC VTLEKCAVCG TACKHLALSD NLKPQEKMFF KSPVETALQY FSHISQVWSF QKKQTDLLIA FYKHRITKLE TAMQEAQQAL VSQDKELSVL RKENGELKKF LAILKESPSR YQGSRSITPR PVGITSPSQS VTPRPSFQHS SQVVSRSSSA ESIPYREAGF GSLGQGGRGL QGRRTPRDSY NETPSPASTH SLSYRTSSAS SGQGIFSRP SPNGHSGHTR VLTPNNFAQR ESTTTLESPL SFQLPVLQTL YQQRHRMGLP SGREAWTTSR
Source	Mammalian cell
Target Names	RNF212B
Protein Names	Recommended name: RING finger protein C14orf164
Expression Region	1-300
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.