



Recombinant *Ashbya gossypii* Elongation factor 1-alpha (TEF)

Product Code	CSB-EP334891DOT-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P41752
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MGKEKTHVNV VVIGHVDSGK STTTGHLIYK CGGIDKRTIE KFEKEAAELG KGSFKYAWVL DKLKAERERG ITIDIALWKF ETPKYHVTVI DAPGHRDFIK NMITGTSQAD CAILIAGGV GEFEAGISKD GQTRHALLA YTLGVKQLIV AINKMDSVKW DESRYQEIVK ETSNFIKKVG YNPKTVPFVP ISGWNGDNMI EATTNAPWYK GWEKETKAGA VKGKTLLEAI DAIEPPVRPT DKALRLPLQD VYKIGGIGTV PVGRVETGVI KPGMVVTFAP SGVTTEVKSVM EMHHEQLEEG VPGDNVGFNV KNVSVKEIRR GNVCGDSKND PPKAAESFNA TVIVLNHPGQ ISAGYSPVLD CHTAHIACKF DELLEKNDRR TGKKLEDSPK FLKAGDAAMV KFVPSKPMCV EAFTDYPPLG RFAVRDMRQT VAVGVIKSVV KSDKAGKVTK AAQKAGKK
Source	<i>E.coli</i>
Target Names	TEF
Protein Names	Recommended name: Elongation factor 1-alpha Short name= EF-1-alpha
Expression Region	1-458
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.