



Recombinant Human Doublecortin domain-containing protein 2C (DCDC2C)

Product Code	CSB-BP006541HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A8MYV0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGTRGPSAPV DTTPAKTIVV YRNGDPFYVG KKFVLSRRRA ATFEALLEQL TEQVDVDFGV RRLFTPTRGH RVLGLDALQA GGKYVAAGRE RFKELDYIHI VPRKPAKIRK LKEIKPVVHC DINVPSKWQT YHRISRHINV FTNGRLFIPP AKIIIPKFSL SDWDIVLATI GEKVFPLGGV RKLFTMNGHL LGDSKDLQDN HFYVAVGLET FKYFPYWKSP RVPSEVQQRY ANVEKNSQRK KKVDSKSGKEP CKYDGIIPKT QDSVYYAKEE KKKTLAEPLV QRGAEGDVYK APTPSKETQG ALDVKEEHNV QLEVPVDQAE SSCDKHCCHN FAYNKHIIKF WLPVILISID FEAVE
Source	Baculovirus
Target Names	DCDC2C
Protein Names	Recommended name: Doublecortin domain-containing protein 2C
Expression Region	1-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.