



Recombinant Human Envoplakin-like protein (EVPLL)

Product Code	CSB-MP007871HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A8MZ36
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<pre>MQASADQVER DILETQKRLQ QDRLNSEQSQ ALQHQQETGS SLKEAEVLLK DLFLDVDKAR RLKHPQAEET EKDIEQLHER VTQECAEYCA LYEKMVLPPR RGIQGR LGTR AGAETEAGLR RPVWAGHGGA GGTDRGAQHR AEGDQRPRRA AAEPGGAGCR HHPEPIPRPT EGGVVARAEP GQPVHALQGC TWQLSALAEQ QRRILQQDWS DLMADPAGVR REYEHFKQHE LLSQEQSVNQ LEEDGKRMVE LRHPAVGPIQ AHQEALKMEW QNFLNLCICQ ETQLQHVEDY SRILCPSSSP H</pre>
Source	Mammalian cell
Target Names	EVPLL
Protein Names	Recommended name: Envoplakin-like protein
Expression Region	1-301
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.