



# Recombinant *Kluyveromyces lactis* SEC14 cytosolic factor (SEC14)

<b>Product Code</b>	CSB-MP335079KKB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P24859
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Kluyveromyces lactis</i> (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) ( <i>Candida sphaerica</i> )
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MVSEQEILES YPQVCPGSL SGTGPNLDSE QEAKLKEFRE LLESGLYKER LDDSTLLRFL RARKFDLEAS KIMYENCEKW RKEFGVDTIF EDFHYEEKPL VAKYYPQYYH KTDNDGRPVY IEELGSVNL T QMYKITTQER MLKNLVWEYE AFVRYRLPAC SRKAGYLVET SCTILDLKGI SISSAAQVLS YVREASNIGQ NYYPERMGKF YLINAPFGFS TAFRLF KPFL DPVTVSKIFI LGSSYQKDLL KQIPAENLPK KFGGQSEVSE AEGGLYLSDI GPWREEEYIG PEGEAPKAFQ L
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SEC14
<b>Protein Names</b>	Recommended name: SEC14 cytosolic factor Alternative name(s): Phosphatidylinositol/phosphatidylcholine transfer protein Short name= PI/PC TP
<b>Expression Region</b>	1-301
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.