



Recombinant *Saccharomyces cerevisiae* DnaJ-related protein SCJ1 (SCJ1)

Product Code	CSB-EP335107SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P25303
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	QDYAILEI DKDATEKEIK SAYRQLSKKY HPDKNAGSEE AHQKFIEVGE AYDVLSDPEK KKIYDQFGAD AVKNGGGGGG PGGPGAGGFH DPFDFIFERMF QGGHGGPGGG FGQRQRQGRP MIKVQEKLSL KQFYSGSSIE FTLNLNDECD ACHGSGSADG KLAQCPDCQG RGVIIQVLRM GIMTQQIQQM CGRCGGTGQI IKNECKTCHG KKVTKKNKFF HVDVPPGAPR NYMDTRVGEA EKGPDFDAGD LVIEFKEKDT ENMGYRRRGD NLYRTEVLSA AEALYGGWQR TIEFLDENKP VKLSRPAHVV VSNGEVEVVK GFGMPKGSKG YGDLYIDYVV VMPKTFKSGQ NMLKDEL
Source	E.coli
Target Names	SCJ1
Protein Names	Recommended name: DnaJ-related protein SCJ1 Short name= J protein SCJ1
Expression Region	22-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.