



Recombinant Human MIF4G domain-containing protein (MIF4GD)

Product Code	CSB-EP013827HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A9UHW6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGEPSREEYK IQSFDAETQQ LLKTALKDPG AVDLEKVANV IVDHSLQDCV FSKEAGRMCY AIIQAESKQA GQSVFRRGLL NRLQQEYQAR EQLRARSLQG WVCYVTFICN IFDYLRVNNM PMMALVNPVY DCLFRLAQP D SLSKEEEVDC LVLQLHRVGE QLEKMNGQRM DELFVLIRDG FLLPTGLSSL AQLLLLEIIE FRAAGWK TTP AAHKYYYSEV SD
Source	E.coli
Target Names	MIF4GD
Protein Names	Recommended name: MIF4G domain-containing protein Alternative name(s): SLBP-interacting protein 1 Short name= hSLIP1
Expression Region	1-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This gene encodes a protein which contains an MIF4G domain.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.