



Recombinant Rice stripe virus Suppressor of RNA silencing p3 (p3)

Product Code	CSB-MP335222REQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P26658
Product Type	Recombinant Protein
Immunogen Species	Rice stripe virus (isolate T) (RSV)
Purity	>85% (SDS-PAGE)
Sequence	MNVFTSSVGS VEFDHPLLE NDLTSLINC DDVHCSSRAL CYIYDIHSSR HPSIDEHQFL RLLHGPDDAV TLGSFLKTLI WILSHDKNLP EEYRLPTIMM SSSYVKFFTE VKPRPPSTNC WTCRMSKDNL PFTVPSVKGF PPDAELYIVP ISDHDGKPKV FDNRKTLYRS PSKKRHKYVI SSDKPPLSAR YVKYVDSSAL ESLPGSSPAV L
Source	Mammalian cell
Target Names	p3
Protein Names	Recommended name: Suppressor of RNA silencing p3 Alternative name(s): Protein NS3 Protein p3
Expression Region	1-211
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.