



Recombinant *Saccharomyces cerevisiae* Protein CWC16 (YJU2)

Product Code	CSB-BP335367SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P28320
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSERKAINKY YPPDYNPLEA EKLSRKMMAK LKTMNKSHAS IRLMTPFSMR CLECNEYIPK SRKFNGKKEL LKEYLDSIK IYRLTISCPR CANSIAFRTD PGNSDYVMEV GGVRNYVPQK PNDDLNAKTA VESIDETLQR LVREKEMEQN EKMGIKEQAD DKMDLLEKRL AKIQQEEDD EELENLRKKK LEMSQRAEMI NRSKHAQQEK AVTTDDLNL VDQVFDNHRQ RTNKPGNNND EKRTPLFNPT STKGKIQKKS SVRTNPLGIV IKRGKSLK
Source	Baculovirus
Target Names	YJU2
Protein Names	Recommended name: Protein CWC16
Expression Region	1-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.