



Recombinant Rice gall dwarf virus Outer capsid protein P8

Product Code	CSB-EP335434RMJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P29077
Product Type	Recombinant Protein
Immunogen Species	Rice gall dwarf virus (RGDV)
Purity	≥85% (SDS-PAGE)
Sequence	MSRQAWIETS ALIECISEYG TKCSFDTFQG LTINDISTLS NLMNQISVAS VGFLNDPRTP LQAMSCEFVN FISTADRHAY MLQKNWFDSD VAPNVTTDNF IATYIKPRFS RTVSDVLRQV NNFALQPMEN PKLISRQLGV LKAYDIPYST PINPMDVARS SANVVGNSVQ RRALSTPLIQ GAQNVTFIVS ESDKIIFGTR SLNPIAPGNF QINVPPWYSD LNVVDARIYF TNSFLGCTIQ NVQVNAVNGN DPVATITVPT DNNPFIVDSV SVVSLSLSGG AINVTAVNL TGYAIAIEGK FNMQMNASPS YYTLSSSLTIQ TSVIDDFGLS AFLEPFRIRL RASGQTEIFS QSMNTLTENL IRQYMPANQA VNIAFVSPWY RFSEARARTIL TFNQPLLPFA SRKLIIRHLW VIMSFIAVFG RYYTVN
Source	E.coli
Protein Names	Recommended name: Outer capsid protein P8 Alternative name(s): Structural protein P8
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.