



Recombinant Simian foamy virus type 1 Protein Bel-2 (bel2)

Product Code	CSB-BP335446SGP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P29170
Product Type	Recombinant Protein
Immunogen Species	Simian foamy virus type 1 (SFVmac) (SFV-1)
Purity	>85% (SDS-PAGE)
Sequence	MNDACLLFPA SKELDSQNVA VILADDMKHS ESILQLGSGA PRKGFPLTSS PLLPMVTPWP FSQDHAAPTL YSLLVAYYKS FQSQKLDPPK WLWQCLGDPS GRKCMVTQFL LPPLGQVRIS CYRNLT SIVI CQAVDPWENN NEADWRKNPM ARPRIKCDHA LCFKVVEGT LWRPHDQKCW LIRLTEGHKY GMEELSPGDW KILQESRPYP YGPIGEDPNL QYAVGVKMKV IGGPLTSTVL ALKALSFHRV NICNMDNPSL GEGHAPLRY S HALKAYGPQY GSCEERVWQT ATKCIGPEEE NYWCEYDHRG FFPMVPNKLS PTWVRHAAPY CIQRFATPYD LQYFANELLP PGFSITTPKG VSYTSDRRLH YGNEGTLQEY NENCDKVKRG YDEISSSDYS DEN
Source	Baculovirus
Target Names	bel2
Protein Names	Recommended name: Protein Bel-2
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.