



Recombinant Rat WAS protein family homolog 1 (Wash1)

Product Code	CSB-MP025972RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B2RYF7
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MTAVKTQHSL AGQVYAVPLI QPDLRREEAI QQVADALQYL QNISGDIFSR ISQRVELSRR QLQAIGERVS LAQAKIEKIK GSKKAIKVFS SAKYPAPEHL QEYNSVFTGA LDPGLQRRPR YRIQSKHRPL DERALQEKLK YFPVCVSTKS EPEDEAEEGL GGLPSNISSI SLLLLFNTTE NLYKKYVFLD PLAGAVTKTH TMLGTEEEKL FDAPLSISKR EQLEQPAPEN YFYVPGLGQV PEIDVPSYLP DLPGVADDLM YSADLGPZIA PSAPGAPEL PAFHTEVAEP FQPEREDGAL LAPPPPPPPP PPPPPAPTA VVSAPQPPMS PDVVTVTGKV AREEDSGSGE AHSASVQGAP KEVVDPSSGR ATLLSIRQA GGIGKAKLRS VKERKLEKKK QKEQEQVRAT SQGGDLMSDL FNKLVMRKKG ISGKGP GTGT SEGPGGAFSR MSDSIPPLPP PQQPAGDEDE DDWES
Source	Mammalian cell
Target Names	Washc1
Protein Names	Recommended name: WAS protein family homolog 1
Expression Region	1-475
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.