



Recombinant Drosophila melanogaster Protein bangles and beads (bnb)

Product Code	CSB-MP335489DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P29746
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MKCQLLFVAT VCLASVWALP VPDEEVAIQP DSGAKAELLT KTVLPTAVEP APLKPEAEKP AETKTIEAKA AVPEQPVDNA NPSPSTPAAE PAKEINSVEL KSSAPDVETA PAIPEKKTLP EEAKPAQENA PVEAEKKQEK TARTEAEPTV EAQPQATKAI EQAPEAPAAAN AEVQKQVVDE VKPQEPKIDA KSAEIPAIPA VVAAEKETPV PEQPARQERI NEIEQKDAKK DAAVAEPAK AAEATPTAAP EAATKSDSNI QVIAPEKCSI ESSPAAAAAS PAAQAAQAKS GEAPKPVDDQ KSTETVAESA PVLKTNAPLA PAGATKVTEA VKEQEKEQPA ADTAAKALPE QKKTEETAAP AGAPEPTAAV APAAVPEAKK IDEAPAAETV VKGEEAIAKP IAQSPSAEPK KSSEKSDKS ESKVDESSES KESEESSESK EN
Source	Mammalian cell
Target Names	bnb
Protein Names	Recommended name: Protein bangles and beads
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.